



SRGAP1 Antibody, HRP conjugated

Product Code	CSB-PA770358LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7Z6B7
Immunogen	Recombinant Human SLIT-ROBO Rho GTPase-activating protein 1 protein (136-192AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	SLIT-ROBO Rho GTPase-activating protein 1 (srGAP1) (Rho GTPase-activating protein 13), SRGAP1, ARHGAP13 KIAA1304
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	SRGAP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.