




# Major pollen allergen Art v 1 Antibody

<b>Product Code</b>	CSB-PA774628HA01AOH
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q84ZX5
<b>Immunogen</b>	Recombinant Artemisia vulgaris Major pollen allergen Art v 1 protein (25-132AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Artemisia vulgaris
<b>Tested Applications</b>	ELISA, WB; Recommended dilution: WB:1:500-1:5000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Major pollen allergen Art v 1 (Defensin-like protein 1) (allergen Art v 1)
<b>Immunogen Species</b>	Artemisia vulgaris
<b>Research Area</b>	Others

## Image



### Western Blot

Positive WB detected in Recombinant protein  
All lanes: Major pollen allergen Art v 1 antibody  
at 3.4µg/ml

### Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution  
predicted band size: 27 kDa  
observed band size: 27 kDa

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.