



TRAPPC11 Antibody, Biotin conjugated

Product Code	CSB-PA801799LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7Z392
Immunogen	Recombinant Human Trafficking protein particle complex subunit 11 protein (704-836AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Trafficking protein particle complex subunit 11, TRAPPC11, C4orf41
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	TRAPPC11
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.