



Ltbp1 Antibody, Biotin conjugated

Product Code	CSB-PA807018LD01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8CG19
Immunogen	Recombinant Mouse Latent-transforming growth factor beta-binding protein 1 protein (20-299AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Latent-transforming growth factor beta-binding protein 1 (LTBP-1) (Transforming growth factor beta-1-binding protein 1) (TGF-beta1-BP-1), Ltbp1
Immunogen Species	Mus musculus (Mouse)
Research Area	Signal Transduction
Target Names	Ltbp1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.