




TRIM59 Antibody

Product Code	CSB-PA808999LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IWR1
Immunogen	Recombinant Human Tripartite motif-containing protein 59 protein (126-181AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tripartite motif-containing protein 59 (RING finger protein 104) (Tumor suppressor TSBF-1), TRIM59, RNF104 TRIM57 TSBF1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	TRIM59

Image



Immunohistochemistry of paraffin-embedded human liver cancer using CSB-PA808999LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.