



TRIM59 Antibody, Biotin conjugated

Product Code	CSB-PA808999LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IWR1
Immunogen	Recombinant Human Tripartite motif-containing protein 59 protein (126-181AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tripartite motif-containing protein 59 (RING finger protein 104) (Tumor suppressor TSBF-1), TRIM59, RNF104 TRIM57 TSBF1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	TRIM59