



TPH2 Antibody, FITC conjugated

Product Code	CSB-PA809000LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IWU9
Immunogen	Recombinant Human Tryptophan 5-hydroxylase 2 protein (1-490AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tryptophan 5-hydroxylase 2 (EC 1.14.16.4) (Neuronal tryptophan hydroxylase) (Tryptophan 5-monoxygenase 2), TPH2, NTPH
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	TPH2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.