



SULF2 Antibody, FITC conjugated

Product Code	CSB-PA810283LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IWU5
Immunogen	Recombinant Human Extracellular sulfatase Sulf-2 protein (463-660AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Extracellular sulfatase Sulf-2 (hSulf-2) (EC 3.1.6.-), SULF2, KIAA1247
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SULF2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.