



LAX1 Antibody, Biotin conjugated

Product Code	CSB-PA811614LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IWW1
Immunogen	Recombinant Human Lymphocyte transmembrane adapter 1 protein (107-298AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Lymphocyte transmembrane adapter 1 (Linker for activation of X cells) (Membrane-associated adapter protein LAX), LAX1, LAX
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	LAX1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.