



MTERF1 Antibody, Biotin conjugated

Product Code	CSB-PA811872LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99551
Immunogen	Recombinant Human Transcription termination factor 1, mitochondrial protein (63-232AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transcription termination factor 1, mitochondrial (Mitochondrial transcription termination factor 1) (mTERF) (mTERF1), MTERF1, MTERF
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	MTERF1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.