




UNC45B Antibody

Product Code	CSB-PA812884LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IWX7
Immunogen	Recombinant Human Protein unc-45 homolog B protein (602-710AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein unc-45 homolog B (Unc-45B) (SMUNC45), UNC45B, CMYA4 UNC45
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	UNC45B

Image



Immunohistochemistry of paraffin-embedded human skeletal muscle tissue using CSB-PA812884LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.