





DYNC2H1 Antibody

Product Code	CSB-PA818763LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NCM8
Immunogen	Recombinant Human Cytoplasmic dynein 2 heavy chain 1 protein (1928-2065AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC, IF; Recommended dilution: WB:1:500-1:5000, IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cytoplasmic dynein 2 heavy chain 1 (Cytoplasmic dynein 2 heavy chain) (Dynein cytoplasmic heavy chain 2) (Dynein heavy chain 11) (hDHC11) (Dynein heavy chain isotype 1B), DYNC2H1, DHC1B DHC2 DNCH2 DYH1B KIAA1997
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	DYNC2H1


Image




Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA818763LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human small intestine tissue using CSB-PA818763LA01HU at dilution of 1:100



Immunofluorescent analysis of HeLa cells using CSB-PA818763LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)



Western Blot

Positive WB detected in Recombinant protein

All lanes: DYNC2H1 antibody at 2.8μg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 45 kDa

Observed band size: 45 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.