



DYNC2H1 Antibody, FITC conjugated

Product Code	CSB-PA818763LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NCM8
Immunogen	Recombinant Human Cytoplasmic dynein 2 heavy chain 1 protein (1928-2065AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cytoplasmic dynein 2 heavy chain 1 (Cytoplasmic dynein 2 heavy chain) (Dynein cytoplasmic heavy chain 2) (Dynein heavy chain 11) (hDHC11) (Dynein heavy chain isotype 1B), DYNC2H1, DHC1B DHC2 DNCH2 DYH1B KIAA1997
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	DYNC2H1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.