



ING5 Antibody, HRP conjugated

Product Code	CSB-PA820185LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WYH8
Immunogen	Recombinant Human Inhibitor of growth protein 5 protein (32-207AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Inhibitor of growth protein 5 (p28ING5), ING5
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	ING5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.