



SAMD8 Antibody, HRP conjugated

Product Code	CSB-PA822254LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96LT4
Immunogen	Recombinant Human Sphingomyelin synthase-related protein 1 protein (1-152AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sphingomyelin synthase-related protein 1 (SMSr) (EC 2.7.8.-) (Ceramide phosphoethanolamine synthase) (CPE synthase) (Sterile alpha motif domain-containing protein 8) (SAM domain-containing protein 8), SAMD8
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SAMD8
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.