



SORCS2 Antibody, FITC conjugated

Product Code	CSB-PA822282LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96PQ0
Immunogen	Recombinant Human VPS10 domain-containing receptor SorCS2 protein (712-786AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	VPS10 domain-containing receptor SorCS2, SORCS2, KIAA1329
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SORCS2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.