



TRPV1 Antibody, HRP conjugated

Product Code	CSB-PA822774LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NER1
Immunogen	Recombinant Human Transient receptor potential cation channel subfamily V member 1 protein (1-155AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transient receptor potential cation channel subfamily V member 1 (TrpV1) (Capsaicin receptor) (Osm-9-like TRP channel 1) (OTRPC1) (Vanilloid receptor 1), TRPV1, VR1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	TRPV1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.