



TRPV1 Antibody, Biotin conjugated

Product Code	CSB-PA822774LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NER1
Immunogen	Recombinant Human Transient receptor potential cation channel subfamily V member 1 protein (1-155AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transient receptor potential cation channel subfamily V member 1 (TrpV1) (Capsaicin receptor) (Osm-9-like TRP channel 1) (OTRPC1) (Vanilloid receptor 1), TRPV1, VR1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	TRPV1