



PF0142 Antibody, HRP conjugated

Product Code	CSB-PA823641LB01FHW
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8U4E6
Immunogen	Recombinant Pyrococcus furiosus Putative L-asparaginase protein (1-175AA)
Raised In	Rabbit
Species Reactivity	Pyrococcus furiosus
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Putative L-asparaginase (EC 3.5.1.1) (L-asparagine amidohydrolase) [Cleaved into: Putative L-asparaginase subunit alpha; Putative L-asparaginase subunit beta]
Immunogen Species	Pyrococcus furiosus (strain ATCC 43587 / DSM 3638 / JCM 8422 / Vc1)
Research Area	Others
Target Names	PF0142
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.