



PF0142 Antibody, FITC conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA823641LC01FHW |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q8U4E6 |
| Immunogen | Recombinant Pyrococcus furiosus Putative L-asparaginase protein (1-175AA) |
| Raised In | Rabbit |
| Species Reactivity | Pyrococcus furiosus |
| Form | Liquid |
| Conjugate | FITC |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Putative L-asparaginase (EC 3.5.1.1) (L-asparagine amidohydrolase) [Cleaved into: Putative L-asparaginase subunit alpha; Putative L-asparaginase subunit beta] |
| Immunogen Species | Pyrococcus furiosus (strain ATCC 43587 / DSM 3638 / JCM 8422 / Vc1) |
| Research Area | Others |
| Target Names | PF0142 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |