



# PF0142 Antibody, Biotin conjugated

<b>Product Code</b>	CSB-PA823641LD01FHW
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q8U4E6
<b>Immunogen</b>	Recombinant Pyrococcus furiosus Putative L-asparaginase protein (1-175AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Pyrococcus furiosus
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	Biotin
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Putative L-asparaginase (EC 3.5.1.1) (L-asparagine amidohydrolase) [Cleaved into: Putative L-asparaginase subunit alpha; Putative L-asparaginase subunit beta]
<b>Immunogen Species</b>	Pyrococcus furiosus (strain ATCC 43587 / DSM 3638 / JCM 8422 / Vc1)
<b>Research Area</b>	Others
<b>Target Names</b>	PF0142
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.