



shha Antibody, Biotin conjugated

Product Code	CSB-PA835612LD01DIL
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92008
Immunogen	Recombinant Zebrafish Sonic hedgehog protein A protein (24-197AA)
Raised In	Rabbit
Species Reactivity	Zebrafish
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sonic hedgehog protein A (SHHA) (Shh unprocessed N-terminal signaling and C-terminal autoprocessing domains) (ShhNC) (VHH-1) [Cleaved into: Sonic hedgehog protein A N-product (Shh N-terminal processed signaling domains) (ShhNp) (Sonic hedgehog protein N-product) (ShhN)], shha, shh vhh1
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Research Area	Others
Target Names	shha
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.