



Ferritin-2, chloroplastic Antibody, FITC conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA836022LC01GGV |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q94IC4 |
| Immunogen | Recombinant Glycine max Ferritin-2, chloroplastic protein (52-257AA) |
| Raised In | Rabbit |
| Species Reactivity | Glycine max |
| Form | Liquid |
| Conjugate | FITC |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Ferritin-2, chloroplastic (EC 1.16.3.1) (SFerH-2) |
| Immunogen Species | Glycine max (Soybean) (Glycine hispida) |
| Research Area | Others |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |