



Ferritin-2, chloroplastic Antibody, Biotin conjugated

Product Code	CSB-PA836022LD01GGV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q94IC4
Immunogen	Recombinant Glycine max Ferritin-2, chloroplastic protein (52-257AA)
Raised In	Rabbit
Species Reactivity	Glycine max
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ferritin-2, chloroplastic (EC 1.16.3.1) (SFerH-2)
Immunogen Species	Glycine max (Soybean) (Glycine hispida)
Research Area	Others
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.