



WDR19 Antibody, HRP conjugated

Product Code	CSB-PA836733LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NEZ3
Immunogen	Recombinant Human WD repeat-containing protein 19 protein (45-226AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	WD repeat-containing protein 19 (Intraflagellar transport 144 homolog), WDR19, IFT144 KIAA1638
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	WDR19
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.