



VTI1A Antibody, HRP conjugated

Product Code	CSB-PA842635LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96AJ9
Immunogen	Recombinant Human Vesicle transport through interaction with t-SNAREs homolog 1A protein (1-175AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Vesicle transport through interaction with t-SNAREs homolog 1A (Vesicle transport v-SNARE protein Vti1-like 2) (Vti1-rp2), VTI1A
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	VTI1A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.