



SLC4A9 Antibody, Biotin conjugated

Product Code	CSB-PA842743LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96Q91
Immunogen	Recombinant Human Anion exchange protein 4 protein (114-402AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Anion exchange protein 4 (AE 4) (Anion exchanger 4) (Sodium bicarbonate cotransporter 5) (Solute carrier family 4 member 9), SLC4A9, AE4 SBC5
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SLC4A9
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.