





ABCA5 Antibody

Product Code	CSB-PA845163LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WWZ7
Immunogen	Recombinant Human ATP-binding cassette sub-family A member 5 protein (1230-1355AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	ATP-binding cassette sub-family A member 5, ABCA5, KIAA1888
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	ABCA5

Image



Immunohistochemistry of paraffin-embedded human liver cancer using CSB-PA845163LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA845163LA01HU at dilution of 1:100