



NEUROD6 Antibody, Biotin conjugated

Product Code	CSB-PA846677LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96NK8
Immunogen	Recombinant Human Neurogenic differentiation factor 6 protein (170-314AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Neurogenic differentiation factor 6 (NeuroD6) (Class A basic helix-loop-helix protein 2) (bHLHa2) (Protein atonal homolog 2), NEUROD6, ATOH2 BHLHA2
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	NEUROD6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.