



# NEUROD6 Antibody, Biotin conjugated

Product Code	CSB-PA846677LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96NK8
Immunogen	Recombinant Human Neurogenic differentiation factor 6 protein (170-314AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Neurogenic differentiation factor 6 (NeuroD6) (Class A basic helix-loop-helix protein 2) (bHLHa2) (Protein atonal homolog 2), NEUROD6, ATOH2 BHLHA2
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	NEUROD6