



LMX1A Antibody, Biotin conjugated

Product Code	CSB-PA848414LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TE12
Immunogen	Recombinant Human LIM homeobox transcription factor 1-alpha protein (9-86AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	LIM homeobox transcription factor 1-alpha (LIM/homeobox protein 1.1) (LMX-1.1) (LIM/homeobox protein LMX1A), LMX1A
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	LMX1A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.