



TIMD4 Antibody, Biotin conjugated

Product Code	CSB-PA850304LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96H15
Immunogen	Recombinant Human T-cell immunoglobulin and mucin domain-containing protein 4 protein (134-301AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	T-cell immunoglobulin and mucin domain-containing protein 4 (TIMD-4) (T-cell immunoglobulin mucin receptor 4) (TIM-4) (T-cell membrane protein 4), TIMD4, TIM4
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	TIMD4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.