



# STK35 Antibody, HRP conjugated

Product Code	CSB-PA851552LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TDR2
Immunogen	Recombinant Human Serine/threonine-protein kinase 35 protein (1-200AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serine/threonine-protein kinase 35 (EC 2.7.11.1) (CLP-36-interacting kinase 1) (CLIK-1) (PDLIM1-interacting kinase 1) (Serine/threonine-protein kinase 35 L1), STK35, CLIK1 PDIK1 STK35L1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	STK35