



STK35 Antibody, Biotin conjugated

Product Code	CSB-PA851552LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TDR2
Immunogen	Recombinant Human Serine/threonine-protein kinase 35 protein (1-200AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serine/threonine-protein kinase 35 (EC 2.7.11.1) (CLP-36-interacting kinase 1) (CLIK-1) (PDLIM1-interacting kinase 1) (Serine/threonine-protein kinase 35 L1), STK35, CLIK1 PDIK1 STK35L1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	STK35
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.