



CEP89 Antibody

Product Code	CSB-PA853505LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96ST8
Immunogen	Recombinant Human Centrosomal protein of 89 kDa protein (25-235AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Centrosomal protein of 89 kDa (Cep89) (Centrosomal protein 123) (Cep123) (Coiled-coil domain-containing protein 123), CEP89, CCDC123
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	CEP89
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.