





UNC80 Antibody

Product Code	CSB-PA854025LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N2C7
Immunogen	Recombinant Human Protein unc-80 homolog protein (124-388AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IF; Recommended dilution: WB:1:500-1:5000, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein unc-80 homolog, UNC80, C2orf21 KIAA1843
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	UNC80

Image



Immunofluorescent analysis of HepG2 cells using CSB-PA854025LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)



Western Blot

Positive WB detected in Recombinant protein

All lanes: UNC80 antibody at 3.5µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 59 kDa

Observed band size: 59 kDa