



UNC80 Antibody, FITC conjugated

Product Code	CSB-PA854025LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N2C7
Immunogen	Recombinant Human Protein unc-80 homolog protein (124-388AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein unc-80 homolog, UNC80, C2orf21 KIAA1843
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	UNC80
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.