



CPT1C Antibody, FITC conjugated

Product Code	CSB-PA855048LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TCG5
Immunogen	Recombinant Human Carnitine O-palmitoyltransferase 1, brain isoform protein (699-803AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Carnitine O-palmitoyltransferase 1, brain isoform (CPT1-B) (EC 2.3.1.21) (CPT IC) (Carnitine O-palmitoyltransferase I, brain isoform) (CPTI-B) (Carnitine palmitoyltransferase 1C), CPT1C, CATL1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	CPT1C
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.