




RB1CC1 Antibody

Product Code	CSB-PA855061LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TDY2
Immunogen	Recombinant Human RB1-inducible coiled-coil protein 1 protein (1544-1594AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IF; Recommended dilution: IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	RB1-inducible coiled-coil protein 1 (FAK family kinase-interacting protein of 200 kDa) (FIP200), RB1CC1, KIAA0203 RBICC
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	RB1CC1

Image



Immunofluorescent analysis of PC-3 cells using CSB-PA855061LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.