



RB1CC1 Antibody, FITC conjugated

Product Code	CSB-PA855061LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TDY2
Immunogen	Recombinant Human RB1-inducible coiled-coil protein 1 protein (1544-1594AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	RB1-inducible coiled-coil protein 1 (FAK family kinase-interacting protein of 200 kDa) (FIP200), RB1CC1, KIAA0203 RBICC
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	RB1CC1