



RB1CC1 Antibody, Biotin conjugated

Product Code	CSB-PA855061LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TDY2
Immunogen	Recombinant Human RB1-inducible coiled-coil protein 1 protein (1544-1594AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	RB1-inducible coiled-coil protein 1 (FAK family kinase-interacting protein of 200 kDa) (FIP200), RB1CC1, KIAA0203 RBICC
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	RB1CC1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.