



NPY6R Antibody, Biotin conjugated

Product Code	CSB-PA857447LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99463
Immunogen	Recombinant Human Putative neuropeptide Y receptor type 6 protein (1-39AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Putative neuropeptide Y receptor type 6 (NPY6-R) (NPY Y1-like receptor) (Putative pancreatic polypeptide receptor 2) (PP2), NPY6R, NPY1RL Y2B
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	NPY6R
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.