



# NPY6R Antibody, Biotin conjugated

Product Code	CSB-PA857447LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99463
Immunogen	Recombinant Human Putative neuropeptide Y receptor type 6 protein (1-39AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Putative neuropeptide Y receptor type 6 (NPY6-R) (NPY Y1-like receptor) (Putative pancreatic polypeptide receptor 2) (PP2), NPY6R, NPY1RL Y2B
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	NPY6R