



PIK3R4 Antibody, Biotin conjugated

Product Code	CSB-PA858716LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99570
Immunogen	Recombinant Human Phosphoinositide 3-kinase regulatory subunit 4 protein (1300-1358AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Phosphoinositide 3-kinase regulatory subunit 4 (PI3-kinase regulatory subunit 4) (EC 2.7.11.1) (PI3-kinase p150 subunit) (Phosphoinositide 3-kinase adaptor protein), PIK3R4, VPS15
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PIK3R4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.