



P2RX4 Antibody, Biotin conjugated

Product Code	CSB-PA859513LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99571
Immunogen	Recombinant Human P2X purinoceptor 4 protein (91-176AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	P2X purinoceptor 4 (P2X4) (ATP receptor) (Purinergic receptor), P2RX4
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	P2RX4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.