





GABARAPL1 Antibody

Product Code	CSB-PA861986LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H0R8
Immunogen	Recombinant Human Gamma-aminobutyric acid receptor-associated protein-like 1 protein (1-117AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Gamma-aminobutyric acid receptor-associated protein-like 1 (Early estrogen-regulated protein) (GABA(A) receptor-associated protein-like 1) (Glandular epithelial cell protein 1) (GEC-1), GABARAPL1, GEC1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	GABARAPL1


Image



Immunohistochemistry of paraffin-embedded human liver cancer using CSB-PA861986LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human brain tissue using CSB-PA861986LA01HU at dilution of 1:100



Immunofluorescent analysis of PC-3 cells using CSB-PA861986LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.