



GNB4 Antibody, Biotin conjugated

Product Code	CSB-PA862056LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9HAV0
Immunogen	Recombinant Human Guanine nucleotide-binding protein subunit beta-4 protein (105-185AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Guanine nucleotide-binding protein subunit beta-4 (Transducin beta chain 4), GNB4
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	GNB4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.