




MALD3 Antibody

Product Code	CSB-PA862991HA01MPS
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9M5X7
Immunogen	Recombinant Apple Non-specific lipid-transfer protein (25-115AA)
Raised In	Rabbit
Species Reactivity	Apple
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Non-specific lipid-transfer protein (LTP) (allergen Mal d 3), MALD3
Immunogen Species	Apple
Research Area	Others
Target Names	MALD3

Image



Western Blot

Positive WB detected in Recombinant protein
All lanes: MALD3 antibody at 3.2µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution
predicted band size: 26 kDa
observed band size: 26 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.