



MALD3 Antibody, Biotin conjugated

Product Code	CSB-PA862991HD01MPS
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9M5X7
Immunogen	Recombinant Apple Non-specific lipid-transfer protein (25-115AA)
Raised In	Rabbit
Species Reactivity	Apple
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Non-specific lipid-transfer protein (LTP) (allergen Mal d 3), MALD3
Immunogen Species	Malus domestica (Apple) (Pyrus malus)
Research Area	Others
Target Names	MALD3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.