



YIPF3 Antibody, FITC conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA863928LC01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q9GZM5 |
| Immunogen | Recombinant Human Protein YIPF3 protein (2-148AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Form | Liquid |
| Conjugate | FITC |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Protein YIPF3 (Killer lineage protein 1) (Natural killer cell-specific antigen KLIP1) (YIP1 family member 3) [Cleaved into: Protein YIPF3, 36 kDa form III], YIPF3, C6orf109 KLIP1 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Cardiovascular |
| Target Names | YIPF3 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |